


1277698



# THE UNITED STATES OF AMERICA

TO ALL TO WHOM THESE PRESENTS SHALL COME:

UNITED STATES DEPARTMENT OF COMMERCE

United States Patent and Trademark Office

*January 25, 2005*

**THIS IS TO CERTIFY THAT ANNEXED HERETO IS A TRUE COPY FROM THE RECORDS OF THE UNITED STATES PATENT AND TRADEMARK OFFICE OF THOSE PAPERS OF THE BELOW IDENTIFIED PATENT APPLICATION THAT MET THE REQUIREMENTS TO BE GRANTED A FILING DATE.**

**APPLICATION NUMBER: 60/556,737**

**FILING DATE: *March 26, 2004***

**RELATED PCT APPLICATION NUMBER: *PCT/US04/43499***



Certified by

Under Secretary of Commerce  
for Intellectual Property  
and Director of the United States  
Patent and Trademark Office

11324 U.S. PTO  
032604

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

**PROVISIONAL APPLICATION FOR PATENT COVER SHEET**

This is a request for filing a PROVISIONAL APPLICATION FOR PATENT under 37 CFR 1.53(c).

Express Mail Label No. EV 389269752 US

INVENTOR(S)					
Given Name (first and middle (if any))		Family Name or Surname		Residence (City and either State or Foreign Country)	
Song-Hua Vijay K.		Ke Mahant		San Diego, CA Murrieta, CA	
Additional inventors are being named on the _____ separately numbered sheets attached hereto					
TITLE OF THE INVENTION (500 characters max)					
Multiplexed Nucleic Acid Analysis With Improved Specificity					
Direct all correspondence to: CORRESPONDENCE ADDRESS					
<input checked="" type="checkbox"/> Customer Number: 34284					
OR					
<input type="checkbox"/> Firm or Individual Name					
Address					
Address					
City		State		Zip	
Country		Telephone		Fax	
ENCLOSED APPLICATION PARTS (check all that apply)					
<input checked="" type="checkbox"/> Specification Number of Pages 9					
<input type="checkbox"/> Drawing(s) Number of Sheets					
<input type="checkbox"/> Application Data Sheet. See 37 CFR 1.76					
<input type="checkbox"/> CD(s), Number					
<input type="checkbox"/> Other (specify)					
METHOD OF PAYMENT OF FILING FEES FOR THIS PROVISIONAL APPLICATION FOR PATENT					
<input checked="" type="checkbox"/> Applicant claims small entity status. See 37 CFR 1.27.					
<input type="checkbox"/> A check or money order is enclosed to cover the filing fees.					
<input checked="" type="checkbox"/> The Director is hereby authorized to charge filing fees or credit any overpayment to Deposit Account Number: 502191					
<input type="checkbox"/> Payment by credit card. Form PTO-2038 is attached.					
FILING FEE Amount (\$) 80.00					
The invention was made by an agency of the United States Government or under a contract with an agency of the United States Government.					
<input checked="" type="checkbox"/> No.					
<input type="checkbox"/> Yes, the name of the U.S. Government agency and the Government contract number are: _____					

19587 U.S. PTO  
60/556737

[Page 1 of 1]

Respectfully submitted,

SIGNATURE



TYPED or PRINTED NAME Martin Fessenmaier

TELEPHONE 714-641-5100

Date 03/26/04

REGISTRATION NO. 46697

(If appropriate)

Docket Number: 100788.0023PRO

**USE ONLY FOR FILING A PROVISIONAL APPLICATION FOR PATENT**

This collection of information is required by 37 CFR 1.51. The information is required to obtain or retain a benefit by the public which is to file (and by the USPTO to process) an application. Confidentiality is governed by 35 U.S.C. 122 and 37 CFR 1.14. This collection is estimated to take 8 hours to complete, including gathering, preparing, and submitting the completed application form to the USPTO. Time will vary depending upon the individual case. Any comments on the amount of time you require to complete this form and/or suggestions for reducing this burden should be sent to the Chief Information Officer, U.S. Patent and Trademark Office, U.S. Department of Commerce, P.O. Box 1450, Alexandria, VA 22313-1450. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Mail Stop Provisional Application, Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450.

If you need assistance in completing the form, call 1-800-PTO-9199 and select option 2.

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

# FEE TRANSMITTAL for FY 2004

Effective 10/01/2003. Patent fees are subject to annual revision.

☒ Applicant claims small entity status. See 37 CFR 1.27

TOTAL AMOUNT OF PAYMENT (\$ ) 80.00

## Complete if Known

Application Number	
Filing Date	March 26, 2004
First Named Inventor	Song-Hua Ke
Examiner Name	
Art Unit	
Attorney Docket No.	100788.0023PRO

## METHOD OF PAYMENT (check all that apply)

☐ Check ☐ Credit card ☐ Money Order ☐ Other ☐ None
☒ Deposit Account:
 Deposit  
Account  
Number  
Deposit  
Account  
Name

502191

Rutan &amp; Tucker

The Director is authorized to: (check all that apply)

☒ Charge fee(s) indicated below ☒ Credit any overpayments☒ Charge any additional fee(s) or any underpayment of fee(s)☐ Charge fee(s) indicated below, except for the filing fee to the above-identified deposit account.

## FEE CALCULATION

## 1. BASIC FILING FEE

Large Entity Fee Code (\$)	Small Entity Fee Code (\$)	Fee Description	Fee Paid
1001 770	2001 385	Utility filing fee	
1002 340	2002 170	Design filing fee	
1003 530	2003 265	Plant filing fee	
1004 770	2004 385	Reissue filing fee	
1005 160	2005 80	Provisional filing fee	80.00

SUBTOTAL (1) (\$ ) 80.00

## 2. EXTRA CLAIM FEES FOR UTILITY AND REISSUE

Total Claims	Extra Claims	Fee from below	Fee Paid
Independent Claims	-20** =	X	
Multiple Dependent	-3** =	X	

Large Entity Fee Code (\$)	Small Entity Fee Code (\$)	Fee Description
1202 18	2202 9	Claims in excess of 20
1201 86	2201 43	Independent claims in excess of 3
1203 290	2203 145	Multiple dependent claim, if not paid
1204 86	2204 43	** Reissue independent claims over original patent
1205 18	2205 9	** Reissue claims in excess of 20 and over original patent

SUBTOTAL (2) (\$ )

\*\*or number previously paid, if greater; For Reissues, see above

## FEE CALCULATION (continued)

## 3. ADDITIONAL FEES

Large Entity Fee Code (\$)	Small Entity Fee Code (\$)	Fee Description	Fee Paid
1051 130	2051 65	Surcharge - late filing fee or oath	
1052 50	2052 25	Surcharge - late provisional filing fee or cover sheet	
1053 130	1053 130	Non-English specification	
1812 2,520	1812 2,520	For filing a request for <i>ex parte</i> reexamination	
1804 920*	1804 920*	Requesting publication of SIR prior to Examiner action	
1805 1,840*	1805 1,840*	Requesting publication of SIR after Examiner action	
1251 110	2251 55	Extension for reply within first month	
1252 420	2252 210	Extension for reply within second month	
1253 950	2253 475	Extension for reply within third month	
1254 1,480	2254 740	Extension for reply within fourth month	
1255 2,010	2255 1,005	Extension for reply within fifth month	
1401 330	2401 165	Notice of Appeal	
1402 330	2402 165	Filing brief in support of an appeal	
1403 290	2403 145	Request for oral hearing	
1451 1,510	1451 1,510	Petition to institute a public use proceeding	
1452 110	2452 55	Petition to revive - unavoidable	
1453 1,330	2453 665	Petition to revive - unintentional	
1501 1,330	2501 665	Utility issue fee (or reissue)	
1502 480	2502 240	Design issue fee	
1503 640	2503 320	Plant issue fee	
1460 130	1460 130	Petitions to the Commissioner	
1807 50	1807 50	Processing fee under 37 CFR 1.17(q)	
1808 180	1808 180	Submission of Information Disclosure Stmt	
8021 40	8021 40	Recording each patent assignment per property (times number of properties)	
1809 770	2809 385	Filing a submission after final rejection (37 CFR 1.129(a))	
1810 770	2810 385	For each additional invention to be examined (37 CFR 1.129(b))	
1801 770	2801 385	Request for Continued Examination (RCE)	
1802 900	1802 900	Request for expedited examination of a design application	

Other fee (specify)

\*Reduced by Basic Filing Fee Paid

SUBTOTAL (3) (\$ )

## SUBMITTED BY

Name (Print/Type)

Martin Fessenmaier

Registration No.  
(Attorney/Agent)

46697

(Complete if applicable)

Telephone 714-641-5100

Signature

Date

March 26, 2004

**WARNING:** Information on this form may become public. Credit card information should not be included on this form. Provide credit card information and authorization on PTO-2038.

This collection of information is required by 37 CFR 1.17 and 1.27. The information is required to obtain or retain a benefit by the public which is to file (and by the USPTO to process) an application. Confidentiality is governed by 35 U.S.C. 122 and 37 CFR 1.14. This collection is estimated to take 12 minutes to complete, including gathering, preparing, and submitting the completed application form to the USPTO. Time will vary depending upon the individual case. Any comments on the amount of time you require to complete this form and/or suggestions for reducing this burden, should be sent to the Chief Information Officer, U.S. Patent and Trademark Office, U.S. Department of Commerce, P.O. Box 1450, Alexandria, VA 22313-1450. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450.

If you need assistance in completing the form, call 1-800-PTO-9199 and select option 2.

# MULTIPLEXED NUCLEIC ACID ANALYSIS WITH IMPROVED SPECIFICITY

## Field of The Invention

Molecular diagnosis and devices therefor.

## 5 Background of The Invention

Human papillomavirus (HPV) is now considered a major cause of cervical cancer, killing more than 200000 women around the world each year. The HPV virus is relatively common and more than 100 distinct types of HPV have been identified, some of which are considered "high-risk" for the development of cancer. Detection of such high-risk types of HPV has significant  
10 impact on diagnosis, prevention, treatment and management of cervical cancer in HPV-infected women.

To date, molecular HPV diagnosis relies on various formats of hybridization technology, including southern blot, dot blot, line blot, and *in situ* hybridization. For example, HybridCapture 2 from Digene is a nucleic acid hybridization microplate assay based on chemiluminescence for  
15 the qualitative detection, and differentiating low-risk from high risk groups. Other commercially available tests employ similar methods and may detect presence of various types of HPV in a patient sample. However, currently known diagnostic methods based on hybridization often lack specificity due to cross-hybridization. Cross-hybridization may result in a false positive signal due to closely related types of HPV (*e.g.*, where a target DNA has only a single or few  
20 mismatches to the probes being used). To overcome problems associated with cross-hybridization, a number of approaches have been developed and used such as stringent hybridization and washing, use of PNA, super G and C, universal base stretch, locked DNA, etc. However, most of such approaches have limitations in the capability to increase specificity and frequently increase cost or automation difficulties.

25 Hybridization specificity may be driven by control of the melting temperature. However,  $T_m$ -specific hybridization will frequently lead to false positive results due to sequence similarity of closely related sequence variations. For example, non-specific signals may be obtained where a first sequence variation increases a  $T_m$  while a second sequence variation decreases the  $T_m$ . In

addition to the above approach, the specificity may be further improved by employing multiple probes (e.g., using a second or third labeled probe for differentiating HPV subtypes.

Therefore, while there are numerous methods for nucleic acid based testing for HPV and other pathogens are known in the art, all or almost all of them suffer from various problems, which are even more aggravated where such analysis is performed in a multiplex environment (e.g., a biochip). Consequently, there is still a need to provide improved methods and compositions for molecular diagnostics.

### **Detailed Description**

The inventors generally contemplate diagnostic methods that include a primer extension step in which hybridization specificity is significantly improved by selective choice of PCR oligos for generation of the amplicon and matching choice of primer extension oligo. It should be recognized that contemplated methods are applicable in a variety of diagnostic tests, and that the nature of the nucleic acid to be tested may vary considerably. Therefore, suitable tests include genotyping of biological specimen, SNP analysis, etc., and may start from DNA as well as from RNA (hn RNA, mRNA, or other). Particularly preferred methods are performed using microarray and/or non-microarray methodologies, and offer higher specificity and sensitivity compared to currently used methodologies.

Typically, most contemplated methods will include a first amplification step in which the target material is amplified in a single (e.g., a single tube PCR) or multiple tubes for a subsequent primer extension analysis step. However, it should be recognized that such first amplification may not be necessary in all cases and it is therefore contemplated that the first amplification may be omitted entirely so long as the primer extension criteria match the test nucleic acid.

### ***Amplicon Generation***

In an especially contemplated aspect of the inventive subject matter, PCR oligos are designed such that two criteria are met: First, the choice of the oligos is dictated by the position of the sequence variation that is to be detected (wherein the sequence variation is characteristic to a particular genotype). Therefore, both primers will flank the sequence variation such that the

amplicon and not the primers will include the sequence variation. Second, the primer sequences must be chosen such that a DNA product is generated in the PCR to which the primer extension oligo can bind at a predefined melting point. Furthermore, where possible, allele specific primer may be employed for generating the extension product.

5

### *Primer Extension*

Again, two conditions must be met by the primer extension oligo: First, the primer extension oligo is selected such that the oligo only binds to the intended target, but not to a similar genotype sequence at the same experimental condition. In most preferred aspects, such selectivity is achieved via  $T_m$  control. Alternatively, or additionally, salinity and/or dielectric constant of the medium may be altered to modulate binding selectivity. Furthermore, agents that  
10 disrupt hydrogen bonding (chaotropic agents) may be added. Second, the primer extension oligo is further selected such that extension will only occur with proper 3'-end base pairing (*i.e.*, only a perfect match at the 3'-end of the extension primer will allow extension).

Therefore, where multiple sequence variation are to be tested, it should be recognized that  
15 the amplicon primers are chosen such that the amplicon provides a hybridization target having a  $T_m$  and sequence that is distinct from the  $T_m$  and sequence of a second hybridization target of a second amplicon. Particularly preferred  $T_m$  differences are at least 1 °C, more typically at least 2 °C, and even more typically at least 4°C, while the sequence difference between first and second extension oligo is at least 1 nucleotide, more preferably at least 2 nucleotides, and most  
20 preferably at least 4 nucleotides.

It is still further preferred that the primer extension oligo will include a unique sequence (preferably with little or no homology [less than 50%] to the amplicon and/or extension product) that is used to hybridize the primer (with or without extension) to a predetermined corresponding nucleic acid that is preferably at a predetermined position on the chip. Addition of such unique  
25 sequence to the 5'-end will allow rapid correlation of a primer or extension signal to a predetermined position on an array. Moreover, it is also preferred that the extension oligo will carry a quantifiable marker (*e.g.*, fluorescent, luminescent, absorbing, or otherwise labeled) that allows normalization of the extension signal.

Primers complying with these conditions have been found to overcome problems with methods relying exclusively on T<sub>m</sub>-specific hybridization, which tend to produce false positive results due to sequence similarity. Contemplated methods were also found to overcome problems associated with methods relying exclusively on correct 3'-match, which often lead to false positive results due to the inability of the extension polymerase to proof proper hybridization conditions towards the 5'-end.

Sometimes, for a specific type of HPV, random sequence variation may occur within the detect region, especially at the region where the 3' terminal of the detection probe binds. This situation will prevent primer extension and cause false negative results. To overcome this phenomenon, a second probe is employed at a different location within the same amplified region. The chance of random sequence variation at two defined location at the same time will be a rare occurrence. In the event the first probe fails to detect the virus, the second probe thus increases the probability of detecting the virus and this will thus minimize false negative results.

Figure 1 below depicts an exemplary flowchart for contemplated methods:





Figure 1


In a typical process according to Figure 1, DNA is isolated using a commercially available kit (*e.g.*, from QIAGEN or Gentra). The typing of HPV DNA is performed using allele specific PCR or primer extension. The subsequent procedure of hybridization of products onto chip is then performed to visualize or otherwise detect the typing results.

5            Figure 2 below provides an example of HPV detection from an actual patient sample. Here, each type of HPV probe is spotted in triplicate with positive and negative controls onto a biofilm chip. The result indicated HPV type 45 infection for the patient and the results were verified by DNA sequencing. Exemplary biofilm chips and devices are described in our  
10            copending applications with the serial numbers PCT/US02/03917, PCT/US01/47991, and 10/346879, which are incorporated by reference herein. Figure 3 (A) illustrates the detection results of HeLa cell infected with HPV 18 and Figure 3 (B) shows the results of Caski cell line infected with HPV 16. As in Figure 2, each type of HPV probe is spotted in triplicate with positive and negative control.






**Figure 3A**




**A. HeLa Cell HPV Detection**

**Figure 3B**



**B. Caski Cell HPV Detection**

Figure 4 depicts the comparison results from the method according to the inventive subject matter (A) and from a conventional method based on hybridization (B) using the same DNA samples isolated from ME180 cells line. Figure 4 (A) indicates the presence of HPV-68 in the ME 180 cell line without ambiguous, but Figure 4 (B) shows positive signal for both HPV-29 and HPV-53 besides HPV-68. The extra erroneous signals for HPV-29 and HPV-53 were the results from cross-hybridization.



**Figure 4**

Thus, especially suitable methods may include the steps of (a) extraction of DNA from patient sample, (b) PCR amplification of a region of HPV genome, (c) allele specific asymmetric PCR or primer extension with primers specific to each type of HPV, and (d) optionally labeling the extension oligo. The preferred polymerase for primer extension is a polymerase capable of recognizing structure and conformational differences at or near the 3' end of a primer-template complex. In a typical detection, hybridization of the labeled primer extension products takes place with capture probes immobilized on a support (the label can be fluorescent, radioactive, chemiluminescent, phosphorescent, electrochemical, or chromogenic).

### **Experimental Procedure**

Step 1: Isolate virus DNA from Pap Smear Samples using QIAamp DNA Blood Mini Kit Handbook from QIAGEN or using Genomic DNA Purification Kit from Gentra.

Step 2: Perform PCR to amplify target sequences in the region of E6 and E7 of HPV genome.

Step 3: Perform HPV Type-Specific Primer Extension Assay with Platinum GenoTYPE Tsp DNA Polymerase from Invitrogen. In this step, detectable fluorescent moieties (such as Cy5-dCTP or Alexa-dCTP) are incorporated into primer extension products.

Step 4: Apply sample onto Chips and perform hybridization with probes immobilized on the surface of the chips (HPV Genotyping Chips from AutoGenomics, Carlsbad, CA)

Step 5: Wash out unincorporated fluorescent moiety and un-hybridized primer extension product and detect signals with Cy3 or Cy5 scanning or other fluorescent detection systems.

### ***Amplicon Oligos***

Based on sequence homology at the DNA level, a set of primers or degenerated primers have been designed. The optimal primers have a T<sub>m</sub> about 52 °C to 64 °C.

### **For HPV-16**

Upstream Primer: 5' GTATATAGAGATGGGAATCC 3'

Downstream Primer: 5' GCCTCTACATAAAACCATCC 3'

For HPV-18

Upstream Primer: 5' GTGTATAGAGACAGTATACC 3'

Downstream Primer: 5' GCTTGTACATAAAACCAGCC 3'

5

*Primer Extension Oligos*

For probes of each HPV type, a unique sequences specific for each type of HPV are designed. These probe can be extended with Platinum GenoTYPE Tsp DNA polymerase (from Invitrogen) or other enzymes which can recognize the 3' end structure, only in the conditions that correct Watson-Crick base pairings are existing at the 3' end of detection probes and templates.

10 The probes can be 8 to 40 bases in length with oligonucleotide tag sequences attached to them. The optimal probes have a Tm about 52 °C to 65 °C, but can be variable from 42 °C to 74 °C.

For HPV-16

Probe Sequence (ASPE): 5' GTTGCAGATCATCAAGAACAC 3'

\*For HPV-18

15

Probe Sequence (ASPE) 5' CGACAGGAACGACTCCAACG 3'

20

It should be noted that contemplated compositions and methods invariably provided conclusive positive or negative results for a given HPV genotype, whereas some of alternative methods solely based on hybridization or primer extension would have provided inconclusive results (typically mistyping or false positive results). Furthermore, most results obtained using contemplated compositions and methods could be independently confirmed by running known assays in the art (which typically required multiple tests in at least some cases to obtain a conclusive result).

Thus, specific embodiments and applications of primer extension analyses with improved specificity have been disclosed. It should be apparent, however, to those skilled in the art that many more modifications besides those already described are possible without departing from the inventive concepts herein. The inventive subject matter, therefore, is not to be restricted except in the spirit of the present disclosure. Moreover, in interpreting the specification, all terms should be interpreted in the broadest possible manner consistent with the context. In particular, the terms "comprises" and "comprising" should be interpreted as referring to elements, components, or steps in a non-exclusive manner, indicating that the referenced elements, components, or steps may be present, or utilized, or combined with other elements, components, or steps that are not expressly referenced.

# Document made available under the Patent Cooperation Treaty (PCT)

International application number: PCT/US04/043499

International filing date: 22 December 2004 (22.12.2004)

Document type: Certified copy of priority document

Document details: Country/Office: US  
Number: 60/556,737  
Filing date: 26 March 2004 (26.03.2004)

Date of receipt at the International Bureau: 03 February 2005 (03.02.2005)

Remark: Priority document submitted or transmitted to the International Bureau in compliance with Rule 17.1(a) or (b)



World Intellectual Property Organization (WIPO) - Geneva, Switzerland  
Organisation Mondiale de la Propriété Intellectuelle (OMPI) - Genève, Suisse